

High affinity immunoglobulin gamma Fc receptor I

Catalog No: #AP83633



Package Size: #AP83633-1 50ug #AP83633-2 100ug #AP83633-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	High affinity immunoglobulin gamma Fc receptor I
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IgG Fc receptor I,Fc-gamma RI,FcRI,CD64,Fcg1
Accession No.	Uniprot ID:P26151
Uniprot	P26151
GeneID	14129;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fcgr1

Note: This product is for in vitro research use only