

Alpha-fetoprotein

Catalog No: #AP83634



Package Size: #AP83634-1 50ug #AP83634-2 100ug #AP83634-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-fetoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Alpha-1-fetoprotein,Alpha-fetoglobulin
Accession No.	Uniprot ID:Q3SZ57
Uniprot	Q3SZ57
GeneID	506011;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AFP

Note: This product is for in vitro research use only