

## Fibroblast growth factor 6

Catalog No: #AP83635



Package Size: #AP83635-1 50ug #AP83635-2 100ug #AP83635-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibroblast growth factor 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FGF-6,Heparin secretory-transforming protein 2,HST-2,Heparin-binding growth factor 6,HBGF-6,Fgf-6,Hstf2
Accession No.	Uniprot ID:P21658
Uniprot	P21658
GeneID	14177;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Fgf6

Note: This product is for in vitro research use only