

Heat shock protein HSP 90-alpha

Catalog No: #AP83637



Package Size: #AP83637-1 50ug #AP83637-2 100ug #AP83637-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock protein HSP 90-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Heat shock 86 kDa,HSP 86,Hsp86,Hspca
Accession No.	Uniprot ID:P82995
Uniprot	P82995
GeneID	103692716;299331;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hsp90aa1

Note: This product is for in vitro research use only