

Glutamate--cysteine ligase catalytic subunit

Catalog No: #AP83640



Package Size: #AP83640-1 50ug #AP83640-2 100ug #AP83640-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate--cysteine ligase catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GCS heavy chain,Gamma-ECS,Gamma-glutamylcysteine synthetase,Gclc
Accession No.	Uniprot ID:P19468
Uniprot	P19468
GeneID	25283;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gclc

Note: This product is for in vitro research use only