

Zinc finger protein GLI2

Catalog No: #AP83644



Package Size: #AP83644-1 50ug #AP83644-2 100ug #AP83644-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc finger protein GLI2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Accession No.	Uniprot ID:P55879
Uniprot	P55879
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GLI2

Note: This product is for in vitro research use only