

Solute carrier family 2, facilitated glucose transporter member 4

Catalog No: #AP83648

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83648-1 50ug #AP83648-2 100ug #AP83648-3 1mg

Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Glucose transporter type 4, insulin-responsive, GLUT-4, GLUT4
Accession No.	Uniprot ID:Q9XT10
Uniprot	Q9XT10
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC2A4

Note: This product is for in vitro research use only