

Insulin-like 3

Catalog No: #AP83650



Package Size: #AP83650-1 50ug #AP83650-2 100ug #AP83650-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Leydig insulin-like peptide, Ley-I-L, Relaxin-like factor, RLF
Accession No.	Uniprot ID:P51461
Uniprot	P51461
GeneID	397024;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: INSL3

Note: This product is for in vitro research use only