

Galanin receptor type 1

Catalog No: #AP83655



Package Size: #AP83655-1 50ug #AP83655-2 100ug #AP83655-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galanin receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GAL1-R,GALNR,GALNR1
Accession No.	Uniprot ID:P47211
Uniprot	P47211
GeneID	2587;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GALR1

Note: This product is for in vitro research use only