

# Gamma-glutamyltranspeptidase 1

Catalog No: #AP83657



Package Size: #AP83657-1 50ug #AP83657-2 100ug #AP83657-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Gamma-glutamyltranspeptidase 1  |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Rt  |
| Other Names        | GGT 1, Gamma-glutamyltransferase 1, Glutathione hydrolase 1, Leukotriene-C4 hydrolase, CD224, Ggt |
| Accession No.      | Uniprot ID:P07314   |
| Uniprot            | P07314  |
| GeneID             | 116568;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ggt1

Note: This product is for in vitro research use only