

Granzyme A

Catalog No: #AP83660



Package Size: #AP83660-1 50ug #AP83660-2 100ug #AP83660-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Granzyme A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Autocrine thymic lymphoma granzyme-like serine protease,CTLA-3,Fragmentin-1,T cell-specific serine protease 1,TSP-1,Ctla-3,Ctla3,Mtsp-1
Accession No.	Uniprot ID:P11032
Uniprot	P11032
GeneID	14938;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gzma

Note: This product is for in vitro research use only