

Glyoxylate reductase/hydroxypyruvate reductase

Catalog No: #AP83661



Package Size: #AP83661-1 50ug #AP83661-2 100ug #AP83661-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glyoxylate reductase/hydroxypyruvate reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GLXR,MSTP035
Accession No.	Uniprot ID:Q9UBQ7
Uniprot	Q9UBQ7
GeneID	9380;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GRHPR

Note: This product is for in vitro research use only