

Eukaryotic translation initiation factor 5A-1

Catalog No: #AP83664



Package Size: #AP83664-1 50ug #AP83664-2 100ug #AP83664-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Eukaryotic translation initiation factor 5A-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	eIF-5A-1,Eukaryotic initiation factor 5A isoform 1,eIF-5A,Rev-binding factor,eIF-4D
Accession No.	Uniprot ID:P63241
Uniprot	P63241
GeneID	1984;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EIF5A

Note: This product is for in vitro research use only