

Glucosylceramidase

Catalog No: #AP83670



Package Size: #AP83670-1 50ug #AP83670-2 100ug #AP83670-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glucosylceramidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Acid beta-glucosidase,Beta-glucocerebrosidase,D-glucosyl-N-acylsphingosine glucohydrolase
Accession No.	Uniprot ID:Q70KH2
Uniprot	Q70KH2
GeneID	449572;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GBA

Note: This product is for in vitro research use only