

Glycogen debranching enzyme

Catalog No: #AP83671



Package Size: #AP83671-1 50ug #AP83671-2 100ug #AP83671-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glycogen debranching enzyme
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Glycogen debrancher
Accession No.	Uniprot ID:Q2PQH8
Uniprot	Q2PQH8
GeneID	479931;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AGL

Note: This product is for in vitro research use only