

Hepatocyte nuclear factor 1-alpha

Catalog No: #AP83673



Package Size: #AP83673-1 50ug #AP83673-2 100ug #AP83673-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatocyte nuclear factor 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HNF-1-alpha,Liver-specific transcription factor LF-B1,LFB1,Transcription factor 1,TCF-1,Hnf-1,Hnf-1a,Tcf1
Accession No.	Uniprot ID:P22361
Uniprot	P22361
GeneID	21405;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hnf1a

Note: This product is for in vitro research use only