

Coagulation factor X

Catalog No: #AP83676



Package Size: #AP83676-1 50ug #AP83676-2 100ug #AP83676-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Coagulation factor X
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Stuart factor, Stuart-Prower factor
Accession No.	Uniprot ID:P00742
Uniprot	P00742
GeneID	2159;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: F10

Note: This product is for in vitro research use only