

# Gremlin-2

Catalog No: #AP83680



Package Size: #AP83680-1 50ug #AP83680-2 100ug #AP83680-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gremlin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cysteine knot superfamily 1, BMP antagonist 2, Protein related to DAN and cerberus, PRDC, Cktsf1b2, Prdc
Accession No.	Uniprot ID:O88273
Uniprot	O88273
GeneID	23893;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Grem2

Note: This product is for in vitro research use only