

Fibroblast growth factor receptor 1

Catalog No: #AP83681



Package Size: #AP83681-1 50ug #AP83681-2 100ug #AP83681-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	FGFR-1,Basic fibroblast growth factor receptor 1,MFR,Proto-oncogene c-Fgr,CD331,Flg
Accession No.	Uniprot ID:Q04589
Uniprot	Q04589
GeneID	79114;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgfr1

Note: This product is for in vitro research use only