

## Filamin-B

Catalog No: #AP83682



Package Size: #AP83682-1 50ug #AP83682-2 100ug #AP83682-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Filamin-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FLN-B,ABP-278,ABP-280 homolog,Actin-binding-like protein,Beta-filamin,Filamin homolog 1,Fh1,Filamin-3,Thyroid autoantigen,Truncated actin-binding protein,Truncated ABP,FLN1L,FLN3,TABP,TAP
Accession No.	Uniprot ID:O75369
Uniprot	O75369
GeneID	2317;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FLNB

Note: This product is for in vitro research use only