

Haptoglobin-related protein

Catalog No: #AP83684



Package Size: #AP83684-1 50ug #AP83684-2 100ug #AP83684-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Haptoglobin-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:P00739
Uniprot	P00739
GeneID	3250;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HPR

Note: This product is for in vitro research use only