

Heat shock protein beta-1

Catalog No: #AP83685



Package Size: #AP83685-1 50ug #AP83685-2 100ug #AP83685-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock protein beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	HspB1,25 kDa IAP,Actin polymerization inhibitor,Heat shock 25 kDa protein,HSP 25,Heat shock 27 kDa protein,HSP 27
Accession No.	Uniprot ID:Q00649
Uniprot	Q00649
GeneID	396227;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPB1

Note: This product is for in vitro research use only