

Glutamate receptor 2

Catalog No: #AP83686



Package Size: #AP83686-1 50ug #AP83686-2 100ug #AP83686-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GluR-2,AMPA-selective glutamate receptor 2,GluR-B,GluR-K2,Glutamate receptor ionotropic, AMPA 2,GluA2,Glur2
Accession No.	Uniprot ID:P19491
Uniprot	P19491
GeneID	29627;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gria2

Note: This product is for in vitro research use only