

Phospholipid hydroperoxide glutathione peroxidase, mitochondrial



Catalog No: #AP83690

Orders: order@signalwayantibody.com

Package Size: #AP83690-1 50ug #AP83690-2 100ug #AP83690-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Phospholipid hydroperoxide glutathione peroxidase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PHGPx, Glutathione peroxidase 4, GPx-4
Accession No.	Uniprot ID:Q9N2J2
Uniprot	Q9N2J2
GeneID	286809;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GPX4

Note: This product is for in vitro research use only