

# Solute carrier family 2, facilitated glucose transporter member 3

Catalog No: #AP83693

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP83693-1 50ug #AP83693-2 100ug #AP83693-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Glucose transporter type 3, brain, GLUT-3, Glut-3, Glut3
Accession No.	Uniprot ID:P32037
Uniprot	P32037
GeneID	20527;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Slc2a3

Note: This product is for in vitro research use only