

Guanylate-binding protein 4

Catalog No: #AP83695



Package Size: #AP83695-1 50ug #AP83695-2 100ug #AP83695-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Guanylate-binding protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GTP-binding protein 3,GBP-3,GTP-binding protein 4,GBP-4,Guanine nucleotide-binding protein 4,Guanylate-binding protein 3,Gbp3
Accession No.	Uniprot ID:Q61107
Uniprot	Q61107
GeneID	55932;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gbp4

Note: This product is for in vitro research use only