

Heparin cofactor 2

Catalog No: #AP83696



Package Size: #AP83696-1 50ug #AP83696-2 100ug #AP83696-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heparin cofactor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Heparin cofactor II,HC-II,Protease inhibitor leuserpin-2,HLS2,Serpin D1,HCF2
Accession No.	Uniprot ID:P05546
Uniprot	P05546
GeneID	3053;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SERPIND1

Note: This product is for in vitro research use only