

# Homeobox protein engrailed-2

Catalog No: #AP83697



Package Size: #AP83697-1 50ug #AP83697-2 100ug #AP83697-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Homeobox protein engrailed-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Gg-En-2,EN-2
Accession No.	Uniprot ID:Q05917
Uniprot	Q05917
GeneID	100859435;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EN2

Note: This product is for in vitro research use only