

## High mobility group protein HMG-I/HMG-Y

Catalog No: #AP83698



Package Size: #AP83698-1 50ug #AP83698-2 100ug #AP83698-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	High mobility group protein HMG-I/HMG-Y
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HMG-I(Y),High mobility group AT-hook protein 1,High mobility group protein A1,High mobility group protein R,HMG1Y
Accession No.	Uniprot ID:P17096
Uniprot	P17096
GeneID	3159;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HMGA1

Note: This product is for in vitro research use only