

High mobility group protein HMG-I/HMG-Y

Catalog No: #AP83698



Package Size: #AP83698-1 50ug #AP83698-2 100ug #AP83698-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                    |   |
|--------------------|---|
| Product Name       | High mobility group protein HMG-I/HMG-Y   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Hu  |
| Other Names        | HMG-I(Y),High mobility group AT-hook protein 1,High mobility group protein A1,High mobility group protein R,HMGIY |
| Accession No.      | Uniprot ID:P17096   |
| Uniprot            | P17096  |
| GeneID             | 3159;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HMGA1

Note: This product is for in vitro research use only