

E3 ubiquitin-protein ligase HACE1

Catalog No: #AP83705



Package Size: #AP83705-1 50ug #AP83705-2 100ug #AP83705-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase HACE1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HECT domain and ankyrin repeat-containing E3 ubiquitin-protein ligase 1,HECT-type E3 ubiquitin transferase HACE1,Kiaa1320
Accession No.	Uniprot ID:Q3U0D9
Uniprot	Q3U0D9
GeneID	209462;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hace1

Note: This product is for in vitro research use only