

## Heat-stable enterotoxin receptor

Catalog No: #AP83709



Package Size: #AP83709-1 50ug #AP83709-2 100ug #AP83709-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Heat-stable enterotoxin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	STA receptor,Guanylyl cyclase C,GC-C,Intestinal guanylate cyclase,GUC2C,STAR
Accession No.	Uniprot ID:P25092
Uniprot	P25092
GeneID	2984;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:GUCY2C

Note: This product is for in vitro research use only