

Heparanase

Catalog No: #AP83712



Package Size: #AP83712-1 50ug #AP83712-2 100ug #AP83712-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heparanase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Endo-gluconidase,Hpa
Accession No.	Uniprot ID:Q6YGZ1
Uniprot	Q6YGZ1
GeneID	15442;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hpse

Note: This product is for in vitro research use only