

Hexokinase-2

Catalog No: #AP83713



Package Size: #AP83713-1 50ug #AP83713-2 100ug #AP83713-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hexokinase-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hexokinase type II, HK II, Muscle form hexokinase
Accession No.	Uniprot ID:P52789
Uniprot	P52789
GeneID	3099;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HK2

Note: This product is for in vitro research use only