

Zinc finger protein GLI2

Catalog No: #AP83716



Package Size: #AP83716-1 50ug #AP83716-2 100ug #AP83716-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc finger protein GLI2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GLI family zinc finger protein 2,Tax helper protein,THP
Accession No.	Uniprot ID:P10070
Uniprot	P10070
GeneID	2736;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GLI2

Note: This product is for in vitro research use only