

Heterogeneous nuclear ribonucleoprotein F

Catalog No: #AP83721



Package Size: #AP83721-1 50ug #AP83721-2 100ug #AP83721-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein F
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	hnRNP F,Hnrpf
Accession No.	Uniprot ID:Q9Z2X1
Uniprot	Q9Z2X1
GeneID	98758;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hnrpf

Note: This product is for in vitro research use only