

Heterogeneous nuclear ribonucleoprotein K

Catalog No: #AP83722



Package Size: #AP83722-1 50ug #AP83722-2 100ug #AP83722-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein K
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	hnRNP K,HNRPK,RCJMB04_24e23
Accession No.	Uniprot ID:Q5ZIQ3
Uniprot	Q5ZIQ3
GeneID	426516;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HNRNPK

Note: This product is for in vitro research use only