

## Gamma-aminobutyric acid receptor subunit alpha-3

Catalog No: #AP83724



Package Size: #AP83724-1 50ug #AP83724-2 100ug #AP83724-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gamma-aminobutyric acid receptor subunit alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GABA(A) receptor subunit alpha-3,Gabra-3
Accession No.	Uniprot ID:P20236
Uniprot	P20236
GeneID	24947;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gabra3

Note: This product is for in vitro research use only