

Fibronectin type III domain-containing protein 3B

Catalog No: #AP83726



Package Size: #AP83726-1 50ug #AP83726-2 100ug #AP83726-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibronectin type III domain-containing protein 3B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Factor for adipocyte differentiation 104,HCV NS5A-binding protein 37,FAD104,NS5ABP37,UNQ2421/PRO4979/PRO34274
Accession No.	Uniprot ID:Q53EP0
Uniprot	Q53EP0
GeneID	64778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FNDC3B

Note: This product is for in vitro research use only