

ICOS ligand

Catalog No: #AP83733



Package Size: #AP83733-1 50ug #AP83733-2 100ug #AP83733-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ICOS ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B7 homolog 2,B7-H2,B7-like protein GI50,B7-related protein 1,B7RP-1,CD275,B7H2,B7RP1,ICOSL,KIAA0653
Accession No.	Uniprot ID:O75144
Uniprot	O75144
GeneID	23308;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ICOSLG

Note: This product is for in vitro research use only