

Guanine nucleotide-binding protein G(s) subunit alpha isoforms short



Catalog No: #AP83736

Orders: order@signalwayantibody.com

Package Size: #AP83736-1 50ug #AP83736-2 100ug #AP83736-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha isoforms short
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alternative gene product encoded by XL-exon,Gnas1
Accession No.	Uniprot ID:Q6R0H6
Uniprot	Q6R0H6
GeneID	14683;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gnas

Note: This product is for in vitro research use only