

Histone H4

Catalog No: #AP83738



Package Size: #AP83738-1 50ug #AP83738-2 100ug #AP83738-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histone H4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:P62806
Uniprot	P62806
GeneID	100041230;319155;319156;319157;319158;319159;319160;319161;320332;326619;326620;69386;97122;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hist1h4a

Note: This product is for in vitro research use only