

Far upstream element-binding protein 1

Catalog No: #AP83743



Package Size: #AP83743-1 50ug #AP83743-2 100ug #AP83743-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Far upstream element-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FBP,D3ErtD330e
Accession No.	Uniprot ID:Q91WJ8
Uniprot	Q91WJ8
GeneID	51886;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fubp1

Note: This product is for in vitro research use only