

Histone deacetylase 2

Catalog No: #AP83746



Package Size: #AP83746-1 50ug #AP83746-2 100ug #AP83746-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone deacetylase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HD2,YY1 transcription factor-binding protein,Yy1bp
Accession No.	Uniprot ID:P70288
Uniprot	P70288
GeneID	15182;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hdac2

Note: This product is for in vitro research use only