

## Interleukin enhancer-binding factor 3

Catalog No: #AP83752



Package Size: #AP83752-1 50ug #AP83752-2 100ug #AP83752-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin enhancer-binding factor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Double-stranded RNA-binding protein 76,DRBP76,M-phase phosphoprotein 4,MPP4,Nuclear factor associated with dsRNA,NFAR,Nuclear factor of activated T-cells 90 kDa,NF-AT-90,Translational control protein 80,TCP80,DRBF,MPHOSPH4,NF90
Accession No.	Uniprot ID:Q12906
Uniprot	Q12906
GeneID	3609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ILF3

Note: This product is for in vitro research use only