

Integrin alpha-L

Catalog No: #AP83754



Package Size: #AP83754-1 50ug #AP83754-2 100ug #AP83754-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-L
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD11 antigen-like family member A,Leukocyte adhesion glycoprotein LFA-1 alpha chain,LFA-1A,Leukocyte function-associated molecule 1 alpha chain,Lymphocyte antigen 15,Ly-15,CD11a,Lfa-1,Ly-15
Accession No.	Uniprot ID:P24063
Uniprot	P24063
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Itgal

Note: This product is for in vitro research use only