

Serine protease hepsin

Catalog No: #AP83759



Package Size: #AP83759-1 50ug #AP83759-2 100ug #AP83759-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine protease hepsin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Transmembrane protease serine 1, TMPRSS1
Accession No.	Uniprot ID:P05981
Uniprot	P05981
GeneID	3249;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HPN

Note: This product is for in vitro research use only