

Hypoxia-inducible factor 3-alpha

Catalog No: #AP83760



Package Size: #AP83760-1 50ug #AP83760-2 100ug #AP83760-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Hypoxia-inducible factor 3-alpha |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | HIF-3-alpha,Basic-helix-loop-helix-PAS protein MOP7,Class E basic helix-loop-helix protein 17,bHLHe17,HIF3-alpha-1,Inhibitory PAS domain protein,IPAS,Member of PAS protein 7,PAS domain-containing protein 7,BHLHE17,MOP7,PASD7 |
| Accession No. | Uniprot ID:Q9Y2N7 |
| Uniprot | Q9Y2N7 |
| GeneID | 64344; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HIF3A

Note: This product is for in vitro research use only