

Hypoxanthine-guanine phosphoribosyltransferase

Catalog No: #AP83763



Package Size: #AP83763-1 50ug #AP83763-2 100ug #AP83763-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hypoxanthine-guanine phosphoribosyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	HGPRT,HGPRTase,HPRT
Accession No.	Uniprot ID:Q6WIT9
Uniprot	Q6WIT9
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HPRT1

Note: This product is for in vitro research use only