

Heat shock factor protein 1

Catalog No: #AP83764



Package Size: #AP83764-1 50ug #AP83764-2 100ug #AP83764-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock factor protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HSF 1,Heat shock transcription factor 1,HSTF 1,HSTF1
Accession No.	Uniprot ID:Q00613
Uniprot	Q00613
GeneID	3297;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSF1

Note: This product is for in vitro research use only